



WIRELESS CAMERA AND DISPLAY RECORDER



- Wireless image transmission
- Display, record and playback functions
- Portable with long battery life
- Compatible with all Airtraq models

Prodol Meditec S.A.
Muelle T Olavarri,
5 Las Arenas
48930 Vizcayo Spain
Phone: +34 944 804690
Email: info@airtraq.com

Airtraq LLC.
800 Sun Park Dr.
Fenton, MO 63026 USA
Phone: +1877-624-7929
Fax: +1636-349-3318
Email: info.usa@airtraq.com

WWW.AIRTRAQ.COM

DISTRIBUTED BY:



Where and When You Need It


Guide channel's anatomical shape delivers the endotracheal tube to the glottic opening

Affordability permits full deployment for **immediate** use in OR, ED, ICU, code carts, pre-hospital (EMS), etc.

Effective for **difficult** and **routine** intubations, e.g., face to face, awake, cervical immobilization, etc.

Disposable to **prevent** cross contamination risks

Easy to learn and **simple** to use



Airtraq Avant

Low cost per use and no capital investment

Quick and easy preparation

Optics fully Isolated from patient

Prevents risk of cross infection



Reusable optics discarded
after 50 intubations



Disposable blade
and eyecup



Ready to use

Airtraq sp

Fully disposable single patient use

Ready to use with no set up

From adults to infants

Specialty models also available for nasal and double lumen



Adult Regular

Adult Small

Pediatric

Infant



Infant Nasal

Naso Tracheal

Double Lumen



Regular: 7.0 to 8.5



Small: 6.0 to 7.5



Docking station recharges battery and
displays remaining service life of optics